

## Cole-Parmer Digital Gear Pump 75211-30



The IP23-rated enclosure protects internal components from damage. Digital optical encoder feedback signal keeps speed control to  $\pm 0.3\%$ .

Rating: Not Rated Yet

**Price**

Sales price \$700.00

[Ask a question about this product](#)

Manufacturer [Cole-Parmer](#)

### Description

The IP23-rated enclosure protects internal components from damage. Digital optical encoder feedback signal keeps speed control to  $\pm 0.3\%$ .

Pump head not included, can be ordered from manufacturer.

### Specifications

Temp range	104°F (40°C) maximum
Duty cycle	Continuous
Motor type	Permanent magnet DC (PMDC), IP23
Power (Amps)	3
Power (Hz)	50/60
Dimensions	7-3/4"W x 7-1/4"H x 11-1/2"D
Motor (hp)	1/10
Motor rpm	60 to 3600
Power (VAC)	115
Brand	Cole-Parmer
Manufacturer number	75211-30

*Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.*

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotequipmentsales.com](mailto:sales@biotequipmentsales.com)