

Cole-Parmer Lab Roller



~16 mm rise/fall plus rotation

Rating: Not Rated Yet

Price

Sales price \$675.00

[Ask a question about this product](#)

Manufacturer [Cole-Parmer](#)

Description

This benchtop roller mixer with parallel rollers offers a combined rocking + rolling motion so that tubes, bottles, or flasks mix with little foaming or aeration. It is used for blood tubes, viscous samples, cell culture, blot washing, and other gentle mixing—not high-speed orbital shaking or magnetic stirring.

Key Features / Specifications

Rocking + rolling motion — ~16 mm rise/fall plus rotation. Max load about 3 kg.

Incubator / cold-room use — Typically 4–60 °C and up to ~80% RH. Continuous quiet operation.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com