

## Corepoint Scientific Model FF201WWW/OMHC Flammable Storage -20C Freezer, from 2023!



20 cu ft capacity

Rating: Not Rated Yet

**Price**

Sales price \$2,750.00

[Ask a question about this product](#)

### Description

The Corepoint Scientific White Diamond Series Flammable Storage Freezers are designed for safe storage of flammable and volatile materials in laboratory settings, adhering to NFPA 45 and 70 and OSHA guidelines for general laboratory work areas. The absence of electrical components inside the unit and a vapor-proof enclosure for the compressor's electrical parts enhance safety by reducing ignition risks. (These units are NOT designed for use in environments where volatile/explosive conditions could potentially exist.)

### Features

One solid swing door, right hinged

Six fixed shelves

Keyed door lock

Magnetic door gasket for positive seal

Leveling legs

Probe access port (0.38")

Powder coat exterior

High density urethane foam cabinet and door insulation

HFC-free refrigerant (R290), foam insulation and packaging

Temp Controller: Microprocessor temperature controller

Digital temperature display with sample simulated probe audible and visual high and low temperature alarms

Remote alarm contacts

#### Specifications

Product Type: Freezer

Exterior Dimensions: 30.75" (L) x 28.5" (W) x 74.63" (H)

Interior Dimensions: 25.88" (L) x 21.75" (W) x 65" (H)

Door Material: Solid

Color: White

Capacity: 20 Cu. Ft

Defrost: Manual

Number of Door: 1

Number of Shelves: 6 (Fixed)

Refrigerant: R290

Display Type: Digital

Adjustable Operating Temperature Range: -15°C – -25°C

Temperature Control Type: Microprocessor

Door Type: Swing

Door Hinged: Right

Door Lock Type: Key

Magnetic Door Gasket: Yes

Remote Alarm: Yes

Leveling Legs: Yes

Probe Access Port: 0.38"

H.P: 1/2 HP

Amps: 1.96 Amps

Hertz: 60 Hz

Voltage: 115V

Certifications: UL/C-UL Listed

Weight: 207 lbs

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)