

Coy Vinyl Anaerobic Airlock Chamber



Flexibility allows deeper, higher reach into chamber.

Rating: Not Rated Yet

Price

Sales price \$13,500.00

[Ask a question about this product](#)

Description

Coy Vinyl Anaerobic chambers are constructed of flexible PVC with padded base, efficient vacuum airlock, and aluminum support frame. The vinyl expands and contracts to allow for changes in internal pressure without expelling the expensive gas mix to room atmosphere. The large oval glove ports allow all users to reach deeper and higher into the chamber for greater storage capacity and ergonomics.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com