

Custom Biogenics Systems CBS Isothermal V-1500 AB Series Cryogenic Nitrogen Storage System



30-Liter LN2 capacity

Rating: Not Rated Yet

Price

Sales price \$6,400.00

[Ask a question about this product](#)

Manufacturer [Custom Biogenic Systems](#)

Description

The V-1500 AB has a liquid nitrogen jacket that allows for storage of samples at liquid nitrogen temperature with no direct liquid nitrogen contact; eliminates the possibility of cross-contamination. Comes with 2301 series controller/monitor and seven 5x5x2" cryo box storage racks. Accommodates all standard rack systems.

Features/Specifications:

Wide lid opening for easy access

30-Liter LN2 capacity

Temperatures down to -190°C

5-Liters/day static evaporation rate, 6 days static holding time

9100 2mL vial maximum storage capacity

Rolling caster wheels

110V

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com