

## Custom Biogenics Systems CBS Isothermal V-1500 AB Series Cryogenic Nitrogen Storage System



30-Liter LN2 capacity

Rating: Not Rated Yet

**Price**

Sales price \$6,400.00

[Ask a question about this product](#)

Manufacturer [Custom Biogenic Systems](#)

### Description

The V-1500 AB has a liquid nitrogen jacket that allows for storage of samples at liquid nitrogen temperature with no direct liquid nitrogen contact; eliminates the possibility of cross-contamination. Comes with 2301 series controller/monitor and seven 5x5x2" cryo box storage racks. Accommodates all standard rack systems.

Features/Specifications:

Wide lid opening for easy access  
30-Liter LN2 capacity  
Temperatures down to -190°C  
5-Liters/day static evaporation rate, 6 days static holding time  
9100 2mL vial maximum storage capacity  
Rolling caster wheels  
110V

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)