

Du Pont Sorvall GSA Rotor



Du Pont Sorvall GSA Rotor 6 places

Rating: Not Rated Yet

Price

Sales price \$1,200.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

Du Pont Sorvall GSA Rotor

Diameter - 31cm (12-1/4")

Mass (weight) - 14.9kg (33lbs)

Angle - 28 degrees

Number of Places - 6

Maximum Angular Velocity (in rev/min) - 13,000/min

Maximum Relative Centrifugal Force (RCF) - 27,600

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com