

## Eberbach Con-Torque Tissue Homogenizer



Eberbach Model 7265 ConTorque Tissue Homogenizer power unit with stand (not shown), 0 - 400 rpm. 115 volts, constant torque at all speeds providing continuous, effortless, positive grinding action. Full torque eliminates the possibility of tissue jamming

Rating: Not Rated Yet

**Price**

Sales price \$800.00

[Ask a question about this product](#)

Manufacturer [Eberbach](#)

### Description

Eberbach Model 7265 ConTorque Tissue Homogenizer power unit with stand (not shown), 0 - 400 rpm. 115 volts, constant torque at all speeds providing continuous, effortless, positive grinding action. Full torque eliminates the possibility of tissue jamming pestle in the grinding vessel. crew chuck through a flexible coupling allowing operation at various angles which reduces breakage of glassware.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**