

Elga PURELAB Ultra Bioscience



This ELGA PURELAB Ultra Bioscience water purification system

Rating: Not Rated Yet

Price

Sales price \$2,250.00

[Ask a question about this product](#)

Manufacturer [Elga](#)

Description

!!!SOLD!!!

This ELGA PURELAB Ultra Bioscience water purification system was manufactured in October of 2004 and guarantees water purity at 18.2MΩ-cm, TOC (RO Pretreatment) is 3-10 ppb. Common applications are molecular biology, microbiology, PCR, electrophoresis, DNA sequencing, mammalian cell culture, cell and culture media preparation, in-vitro fertilization, & monoclonal antibody production.

Please call us at (650) 871-5707 to check availability of this product!