

## Eppendorf 5425 Centrifuge





24-place microcentrifuge

Rating: Not Rated Yet

**Price**

Sales price \$1,975.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

This 24-place microcentrifuge offers a broad range of molecular and cell biology applications in Eppendorf and PCR tubes from 0.2 to 5 mL.

Adjustable acceleration and braking ramps; 3 program buttons; short function without continuous pressing of the key. OptiBowl design for quiet operation, even without a rotor lid. Soft-touch one-finger lid closure for ergonomic operation. Aerosol-tight Eppendorf QuickLock lids for simple and safe handling. Noise level is <51 dBA

### Specifications

Display: Large, Brightly Lit LCD

Electrical Requirements: 120 V, 50/60 Hz (US)

Depth: 15.4 in. / 39 cm

With Rotor FA-24x2

Rotary knobs.

Max speed (RCF): 21,300 × g

Max RPM: 15,060

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)