

## Eppendorf 5430



Delux microcentrifuge can accept 8 different rotors, is microprocessor controlled w/up to 50 stored programs, NEW - order rotor from manufacturer

Rating: Not Rated Yet

**Price**

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

**!!!SOLD!!!**

This brand NEW Eppendorf 5430 microcentrifuge offers versatility by accomodating eight different rotors to spin 1.5/2.0mL tubes, 13/16/26/29mm round-bottom tubes, 15/50mL conical tubes, PCR† tube strips, and microplates. They feature menu-driven operation with large backlit LCD, memory for 50 nameable programs, five keys for quick access to routine programs, automatic recognition of rotor and speed limitation, and automatic detection of rotor imbalance. Functions include SOFT-brake for protection of delicate samples, speed programmable short spin for protection of fragile tubes or delicate samples, and end-of-run alarm with six settings for volume. RCF and speed are adjustable in increments of 100. Centrifuges are coldroom compatible.

**Please call us at (650) 871-5707 to check availability of this product!**