

Eppendorf 5810 Benchtop Centrifuge with S-4-104 Rotor



14,000 rpm centrifuge maximum speed, 20,800xg centrifuge maximum g-force

Rating: Not Rated Yet

Price

Sales price \$5,950.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

14,000 rpm centrifuge maximum speed, 20,800xg centrifuge maximum g-force, 4x750mL maximum capacity with S-4-104 rotor, features microprocessor control, adjustable speed short-spin, 10 acceleration and braking rates, 1-99 minute timer or hold, comes with Eppendorf S-4-104 swinging bucket rotor with x4 round buckets and x4 500mL conical bottle adapters, 3,900 rpm maximum speed with S-4-104 rotor, 120V, 50/60Hz, 950W.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com