


Eppendorf 5810 Centrifuge



With A-4-62 rotor.

Rating: Not Rated Yet

Price

Sales price \$4,350.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Benchtop centrifuge with A-4-62 Rotor (swing-bucket rotor with 4x250-ml rectangular buckets)

Max. Rotational Speed: 4,000 rpm

Max. Centrifugal Force: 3,220 x g

Max Load per Bucket: 620 g (includes adapter, tube and contents)

Dimensions (W x D x H): 700 mm x 608 mm x 345 mm (27.6" x 23.9" x 13.6")

Run time: 1 min to 99 min, adjustable in 1 min increments

Max. density of the material to be centrifuged: 1.2 g-ml⁻¹

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com