

## Eppendorf 5810R Centrifuge, 15 amp version



## Centrifuges: Eppendorf 5810R Centrifuge, 15 amp version

---



Slightly worn bowl; with A-4-62 Rotor.

Rating: Not Rated Yet

**Price**

Sales price \$3,650.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Benchtop refrigerating centrifuge with A-4-62 Rotor (swing-bucket rotor with 4x250-ml rectangular buckets). Slightly worn bowl. 15 amp version.

Max. Rotational Speed: 4,000 rpm

Max. Centrifugal Force: 3,220 x g

Max Load per Bucket: 620 g (includes adapter, tube and contents)

Dimensions (W x D x H): 700 mm x 608 mm x 345 mm (27.6" x 23.9" x 13.6")

Run time: 1 min to 99 min, adjustable in 1 min increments

Max. density of the material to be centrifuged: 1.2 g-ml<sup>-1</sup>

Temperature range: -9 °C to 40 °C

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)