

Eppendorf 5910Ri Refrigerated Benchtop Centrifuge w/ S-4xUniversal Rotor, Adapters and FA-6x50 Rotor



14,000 rpm; 22132xg; 4x1000mL; -11°C to 40°C

Rating: Not Rated Yet

Price

Sales price \$14,450.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

The 5910Ri is a highly versatile multipurpose centrifuge with an intuitive touch screen interface and high capacity for a broad range of consumables. The 7" VizioNize touch screen interface provides a documentation function and unique user management feature supporting quality control with multiple authorization levels, password protection, and documentation of all runs and user activities.

Specifications/Features:

- 14,000 rpm maximum speed
- 22132xg maximum g-force
- 4x1000mL maximum capacity,
- temperature control from -11°C to 40°C
- 10/10 acceleration and braking rates
- 10s-99hr59min timer or hold
- microprocessor control and touch screen VizioNize operating interface
- fast temp function
- automatic rotor recognition with rotational speed limit
- automatic rotor imbalance detection
- standby refrigeration keeps rotor at set temperature while lid is closed
- soft touch lid enables easier lid closure

double gas spring lid support
120V, 50/60Hz, 12A

Included:

Eppendorf S-4xUniversal Rotor with 4 buckets (4,500 rpm maximum speed with S-4xUniversal rotor)
(x4) 15mL conical tube adapters
(x4) 50mL conical tube adapters
(x4) 1-Liter bottle adapters
Eppendorf FA-6x50 fixed angle rotor for 50mL conical tubes (12,100 rpm maximum speed with FA-6x50 rotor)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com