


## Eppendorf 5910Ri Refrigerated Benchtop Centrifuge w/ S-4xUniversal Rotor



14,000 rpm maximum speed, 22132xg maximum g-force, 4x1000mL maximum capacity

Rating: Not Rated Yet

**Price**

Sales price \$12,350.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Model: Eppendorf 5910Ri

14,000 rpm maximum speed, 22132xg maximum g-force, 4x1000mL maximum capacity.

Features microprocessor control and touchscreen operating interface, automatic rotor recognition with rotational speed limit, automatic rotor imbalance detection, 10/10 acceleration and braking rates, 10s-99hr59min timer or hold, temperature control from -11°C to 40°C.

With fast temp function. Standby refrigeration keeps rotor at set temperature while lid is closed. Soft touch lid enables easier lid closure. Double gas spring lid support. Comes with Eppendorf S-4xUniversal Rotor with 4 buckets, 4,500 rpm maximum speed with S-4xUniversal Rotor, 120V, 50/60Hz, 12A.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)