

Eppendorf F-45-18-11



Eppendorf F45-18-11-Kit with Lid, 18 place, 1.5 or 2.0mL tubes

Rating: Not Rated Yet

Price

Sales price \$350.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Eppendorf Kit rotor® (F – 45 – 18 – 11 – Kit) for 18 x spin columns and 1.5/2.0 mL tubes.

For use with Eppendorf 5424 and 5424R centrifuges.

Eppendorf Kit rotor® for 18 x spin columns and 1.5/2.0 mL tubes with speed up to 18,111 x g/15,000 rpm.

This rotor is designed with an extended rim that protects tube caps from shearing off.

Additional:

This unique rotor offers a solution for a common problem when working with MiniPrep spin column kits: open tube lids that shear off during the elution step. This problem traditionally occurs because standard rotors do not provide enough support for open tube lids during centrifugation. The Eppendorf Kit rotor is designed with an extended rim that serves as a rest for tube caps during centrifugation and protects them from shearing off. Your tube lids will not tear off anymore.

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com