


Eppendorf F45-24-11 Fixed Angle Rotor



Eppendorf F45-24-11 Fixed-Angle Rotor

Rating: Not Rated Yet

Price

Sales price \$325.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Eppendorf F45-24-11 Fixed-Angle Rotor for use with Eppendorf Centrifuges 5702, 5702R, 5702RH, 5415D, 5415R.

Specifications -

- .Capacity: 24 x 1.5/2.0 mL tubes
- .Max. Rpm: 13,200 rpm
- .Max. rcf: 16,110 x g
- .Autoclavable (121°C, 20 minutes)

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com