

Eppendorf Mini Spin plus



Eppendorf 5453 MiniSpin Plus Centrifuge with Rotor F45-12-11, 120V., 50/60 Hz

Rating: Not Rated Yet

Price

Sales price \$700.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Eppendorf MiniSpin® plus -

Max. speed of 14,500 rpm; max. rcf 14,100 x g

Timer can be set up to 99 minutes

Continuous centrifugation possible

Exceptionally small footprint

Extremely quiet

Low noise levels

Maintenance-free drive

Easy-to-follow digital display

Low sample heating (only 12°C after 20 min. at max. speed)

Separate short-spin key

Centrifuges: Eppendorf Mini Spin plus

Automatic lid release

Built-in power unit

Capacity: 12 x 1.5/2 ml tubes

Autoclavable rotor (121°C, 20 min.)

Acceleration and braking times

Stainless steel rotor lid