

## Eppendorf New Brunswick Excella E25R Large Capacity Top-Loading Refrigerating Incubator Shaker with Universal Platform - Lightly Used from 2022



High-capacity chamber for flasks up to 6 L

Rating: Not Rated Yet

### Price

Sales price \$6,500.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Eppendorf New Brunswick Excella E25R Refrigerated Incubator Shaker

Large-capacity, top-loading refrigerated incubator shaker with universal platform, designed for reliable temperature-controlled shaking applications in microbiology, cell culture, protein expression, and bioprocess laboratories. This lightly used 2022 unit features the trusted New Brunswick orbital drive system, refrigerated temperature control, and high-capacity chamber for flasks up to 6 L. Includes universal platform as shown in photos.

### Key Features

Large Capacity Platform – Accommodates multiple flasks and vessels up to 6 L capacity with included universal platform.

Refrigerated Temperature Control – Operates from 15°C below ambient (minimum 4°C) up to 60°C for versatile incubation applications.

Variable Shaking Speed – Adjustable 50–400 rpm with 1-inch orbital motion for consistent aeration and mixing performance.

Top-Loading Design with Viewing Window – Ergonomic lift-top lid and large viewing window provide easy sample access and monitoring.

Precision Microprocessor Controls – Digital control system with programmable timer, alarms, smooth acceleration/deceleration, and automatic compensation for load and temperature changes.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)