

## Eppendorf / New Brunswick Innova 42R Refrigerated Incubator Shaker



Shaking speed range from 25-400 rpm, temperature control from 20°C below ambient to 80°C, 0.1-99.9 hour timer

Rating: Not Rated Yet

**Price**

Sales price \$7,500.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Mfg no. M1335-0014, Shaking speed range from 25-400 rpm, temperature control from 20°C below ambient to 80°C, 0.1-99.9 hour timer, features triple-eccentric counterbalanced drive, solid state DC brushless motor, 4 storable programs with 15 programmable steps each, 1in. orbital throw, 18x18in. shaker platform dimensions, comes with optional culture drawer, 120V, 60Hz, 12.5A.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**[Biotech Equipment Sales, Inc.](#)**

226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)