

Eppendorf Research Plus 100 12 Channel Pipettor



Current model.

Rating: Not Rated Yet

Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

The Eppendorf Research plus 100 is a variable-volume, 12-channel mechanical air-cushion pipette for simultaneous transfer into 96-well plates (9 mm tip spacing), with a yellow control button, spring-loaded tip cones, and removable individual channels.

Volume range: 10–100 μ L.

Compatible with universal tips.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com