

Eppendorf Research Plus 1000 Single Channel Pipettor



Current model.

Rating: Not Rated Yet

Price

Sales price \$150.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

The Eppendorf Research plus 1000 is a variable-volume, single-channel mechanical air-cushion pipette for accurate, ergonomic manual pipetting of aqueous solutions, with a blue control button and a spring-loaded tip cone.
Volume range: 100–1000 μ L.

Compatible with universal tips.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com