

Eppendorf Research Plus 20 Single Channel Pipettor



Current model.

Rating: Not Rated Yet

Price

Sales price \$150.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

The Eppendorf Research plus 20 is a variable-volume, single-channel mechanical air-cushion pipette for precise, ergonomic manual pipetting of small aqueous volumes, with a spring-loaded tip cone and either a light-gray or yellow control button depending on the tip system. Compatible with universal tips.

Volume range: 2–20 μ L.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com