

Eppendorf Research Plus 5mL Single Channel Pipettor



Current model.

Rating: Not Rated Yet

Price

Sales price \$225.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

The Eppendorf Research plus 5 mL is a variable-volume, single-channel mechanical air-cushion pipette for larger aqueous volumes, with a lilac control button, a fully autoclavable design, and a secondary adjustment for non-aqueous liquids or altitude.

Volume range: 0.5–5 mL.

It uses epT.I.P.S. in the 0.1–5 mL range (about 120 mm long).

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com