

Eppendorf Swing-Bucket Rotor s-4-72



In box. 3,900 RPM; RCF of 3,214 x g

Rating: Not Rated Yet

Price

Sales price \$2,450.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Designed for use with the Eppendorf 5810/R Series centrifuges. Maximum speed of 3,900 RPM and a maximum relative centrifugal force (RCF) of 3,214 x g. It can accommodate a maximum capacity of up to 4 x 750ml, making it suitable for a range of applications in biological, biochemical, and clinical research. A reliable and versatile centrifugation accessory ideal for high-throughput sample processing in laboratories.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

