

## Eppendorf ThermoMixer F1.5 with 24-place head and ThermoTop



With thermoblock for 24 reaction vessels, 1.5 mL.

Rating: Not Rated Yet

### Price

Sales price \$1,975.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Heating and mixing in 1.5 mL reaction vessels. Efficient mixing to 1,500 rpm. Outstanding 2DMix-Control technology results in excellent mixing performance. Anti-spill technology prevents lid wetting and cross-contamination. Compatible with Eppendorf ThermoTop (condens.protect®) for condensation-free sample handling. Predefined temperature keys (37 °C, 42 °C, 56 °C, and 95 °C) for simple and intuitive operation. Very fast heat up speed of up to 11 °C/min, reducing your waiting time.

100 – 130 V/50 – 60 Hz (US/JP/South America/TW)

Specifications Eppendorf ThermoMixer® F1.5

Product type: ThermoMixer

Interfaces: USB interface

Power supply: 100 – 130 V, 50 – 60 Hz

Max. power consumption: 200 W

Dimensions (W × D × H): 20.6 × 30.4 × 17.0 cm / 8.1 × 11.9 × 6.7 in

Product weight: 6.3 kg / 13.9 lb

Max. speed: 1,500 rpm

Temperature range: Min: 4 °C above RT, max: 100 °C

Mixing frequency: 300 – 1,500 rpm

Mixing and vortexing radius: 1.5 mm

Temperature settings: 1 °C / 100 °C

Mixing orbit: 3 mm

Temperature accuracy: Max. ±0.5 °C at 20 – 45 °C

Heating rateMax. 11 °C/min

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)