

## Fischer Scientific AccuSpin Micro 17



Capacity: 48mL (24 x 1.5/2.0mL tubes) Maximum speed/RCF: 13,300rpm/17,000xG Ventilated or refrigerated models Easy-to-use controls with digital display Adjustable 99-min. timer and quick-spin feature for short runs Memory retains last setting for rep

Rating: Not Rated Yet

**Price**

Sales price \$1,100.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

- Capacity: 48mL (24 x 1.5/2.0mL tubes)
- Maximum speed/RCF: 13,300rpm/17,000xG
- Ventilated or refrigerated models
- Easy-to-use controls with digital display
- Adjustable 99-min. timer and quick-spin feature for short runs
- Memory retains last setting for repeatable runs