

## Fisher FBG45RSLA 2-Sliding-Glass Door Digital Refrigerator from 2020



45 cu. ft. (1274L) upright laboratory refrigerator

Rating: Not Rated Yet

**Price**

Sales price \$3,400.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

The Fisher Scientific FBG45RSLA is a 45 cu. ft. (1274L) upright laboratory refrigerator designed for semi-critical sample storage in non-corrosive, non-flammable environments. Featuring a vapor compression system with eco-friendly R290 refrigerant and adaptive defrost, it maintains a temperature range of 2°C to 10°C with a 5°C setpoint. The unit includes eight adjustable wire shelves, a digital eye-level controller, and safety features like remote alarm terminals and a power failure alarm. Certified by cULus and ENERGY STAR, it operates quietly at 47.5 dBA, making it ideal for indoor lab use.

### Specifications:

Capacity: 45 cu. ft. (1274L)

Temperature Range: 2°C to 10°C (5°C setpoint)

Refrigerant: R290 (85g, flammable)

Defrost: Adaptive

Shelves: 8 adjustable wire shelves, 45 kg (100 lbs.) capacity each

Dimensions (H x W x D): Interior: 61.75 x 49 x 27 in. (1569 x 1245 x 686 mm); Exterior: 81 x 52.7 x 32 in. (2057 x 1339 x 813 mm)

Electrical Requirements: 115V, 60Hz, 5A, NEMA 5-15P plug, 6 ft. cord

Sound Level: 47.5 dBA  
Weight: ~255 kg (561 lbs.)  
Certifications: cULus, ENERGY STAR, Prop 65

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)