

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Fisher GPD 15D Dual Water Bath with 10L and 5L Baths



5L and 10L capacities with separate digital controls

Rating: Not Rated Yet

Price

Sales price \$895.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

The Fisher Scientific Isotemp GPD 15D (model FSGPD15D) is a deluxe dual-chamber general-purpose digital water bath designed for laboratory use. It features two independent stainless steel chambers — a 5L (smaller) and a 10L (larger) — allowing simultaneous operation at different temperatures. This compact, reliable unit provides precise temperature control from ambient to 100°C, making it ideal for incubating samples, enzyme reactions, cell culture, microbiology, and other routine lab applications.

This unit offers excellent value for labs needing flexible, space-efficient incubation capacity.

Key Features and Specifications

Dual Independent Chambers — 5L and 10L capacities with separate digital controls, enabling two different temperatures/protocols to run at the same time

Temperature Range & Precision — Ambient to 100°C (with lid); stability $\pm 0.1^\circ\text{C}$, uniformity $\pm 0.2^\circ\text{C}$ at 37°C, and 0.1°C resolution/display.

Includes 3 user-defined temperature presets

Digital Controller with Advanced Features — Easy-to-use interface with timer, audible alarms, auto on/off, Quik-CAL calibration, and over-temperature safety protection

Stainless Steel Construction — Corrosion-resistant 304 stainless steel tanks with rounded corners for easy cleaning; includes hinged see-through or stainless lids

Compact & Practical Design — Overall dimensions approx. 15.4" x 23.1" x 9.2" (39.2 x 58.7 x 23.3 cm); dual heaters (300W + 800W); 100-115V or 200-230V, 50/60 Hz operation. Drain ports and safety features included.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc. has been based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com