

## Fisher Isotemp 1016S



Fisher Isotemp 1016S 6L upright refrigerating heating recirculating bath, digital set/read

Rating: Not Rated Yet

**Price**

Sales price \$950.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

**!!!SOLD!!!**

This Fisher Isotemp 1016S 6L upright refrigerating heating recirculating bath has a digital set/read control, -20 to 150C temp range, 210W of cooling at 20C, and is equipped overtemp & low-level cutouts with visual alarm.

**Please call us at (650) 871-5707 to check availability of this product!**