

Fisher Science Education Portable Balance



Maximum capacity: 220 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

The Fisher Science Education classroom scale is a 220 g compact model, an economical portable scale for general classroom weighing. It has a slim stackable ABS body, a large LCD, two-button operation, keypad calibration with an external weight.

Key specifications

Maximum capacity: 220 g

Readability: 0.1 g

Stabilization time: 2 s (tare to capacity by subtraction; pan about 125 × 132 mm)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com