


## Fisher Scientific GDP 15 D (Model FSGPD15D) Water Bath



5L & 10L Dual Deluxe General Purpose Water Bath

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

Model: FSGPD15D

This 5L & 10L Dual Deluxe General Purpose Water Bath is a compact water bath for industrial, clinical pharmaceutical or biomedical applications.

### Features:

- Smaller footprint, compared to previous models
- Advanced temperature controller is designed for extended functionality
- New audible alarms
- Conveniently save commonly used settings with temperature presets
- Simple to operate and monitor with icon-based graphical display
- Enhanced hinge design supports open lid configuration and easy lid removal
- Help prevent bath damage and overheating with low-fluid protection
- Compatible with Thermal beads (sold separately)

### Specifications:

- Chamber Material: Stainless Steel
- Description: 5L & 10L Dual Deluxe General Purpose Water Bath
- Display: Monochrome LCD
- Includes: Polycarbonate Gable Cover, Diffuser Tray
- Model: Isotemp GPD 15D
- Temperature Stability:  $\pm 0.2^{\circ}\text{C}$
- Dimensions (D x W x H) Interior: 6.1 x 11.8 x 5.99 in. (154 x 300 x 150 mm); 11.9 x 13 x 5.99 in. (301 x 330 x 150 mm)
- Amperage: 8.2-9.7 / 3.8-4.7 (120V/230V)

Controller Type: Digital  
Capacity (Metric): 15 L (5 L/10 L)  
Temperature Range: Ambient to 100°C  
Heating Capacity: 1100 w  
Material: Epoxy Powder-coated Steel  
Temperature Presets: 3  
Temperature Uniformity:  $\pm 0.2^{\circ}\text{C}$   
Dimensions (L x W x H): 15.4 x 23.1 x 9.2 in. (392 x 587 x 233 mm)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)