

Fisher Scientific MaximaDry



Fisher Scientific MaximaDry corrosion resistant oil-free vacuum pump

Rating: Not Rated Yet

Price

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

!!!SOLD!!!

This Fisher Scientific MaximaDry corrosion resistant oil-free vacuum pump, ideal for use with rotary evaporators, gel dryers, filtration flasks and manifolds, vacuum ovens and concentrators. Quiet, reliable operation without the hazards of handling pump oil or contaminated water. Excellent replacement for rotary vane pumps and water aspirators, especially in applications where corrosive vapors are encountered. 2-Stage; 115V, 60Hz; 1/3hp; 15psig; 35L/min.; 6 torr. Manufactured in 2007.

Please call us at (650) 871-5707 to check availability of this product!