

## Fisher Scientific Sonic Dismembrator Model 150 E



Maximum Output: 150 Watts w/fine tip generator

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

Model: 150E

This Fisher Scientific 150 Series Digital Sonic Dismembrator offers continuous or pulsed operation, up to 150 Watts at the horn tip, and a frequency of 40 kHz. There is also an LED display that helps control parameter settings and continuous monitoring.

### Specifications:

Maximum Output: 150 Watts  
Input Power: 117 V, 50/60 Hz, 2.0A  
Height: Approximately 7.9"  
Width: Approximately 13.4"  
Depth: Approximately 9.0"

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

