


Fisher Scientific UV Crosslinker



Five 254nm Shortwave Bulbs; two modes of operation - Energy or Time

Rating: Not Rated Yet

Price

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

Cross linking DNA or RNA to membranes, in preparation for multiple probings of southern's, northern's, dot or slot blots.

The unit comes equipped with five 254nm Shortwave Bulbs
Two modes of operation - Energy or Time

Frequency: 50/60 Hz; Depth (English)14.37 in.; 115 V

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com