

Fisherbrand Basic 60 L Gravity Incubator



Temperature — Ambient +7 °C to 75 °C.

Rating: Not Rated Yet

Price

Sales price \$1,100.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

Fisherbrand Basic 60 L Gravity Incubator is a microbiological incubator. Gravity convection (no fan) gives gentle heat for cultures and warmed storage.

Key Features / Specifications

Chamber — 60 L. Stainless-steel interior. Interior about 16.3 × 13.9 × 20.0 in (414 × 354 × 508 mm). Exterior about 22.2 × 20.9 × 28.3 in (565 × 530 × 720 mm). Weight ~40 kg / 88 lb. Stackable without a kit.

Temperature — Ambient +7 °C to 75 °C. Single-setpoint digital control (0.1 °C steps). At 37 °C: stability about ±0.3 °C, spatial uniformity about ±0.8–0.9 °C. Over-temperature cutoff.

Doors and access — Solid outer door plus inner glass door so you can check samples with less heat loss. Standard access port. Gravity (natural) convection.

Shelves — as pictured.

Electrical / certification — Matches the plate: 120 V, 60 Hz, 2.5 A, 300 W, NEMA 5-15. cCSAus. Allow rear clearance (~80 mm). Not a CO₂ incubator.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com