

Fisherbrand Isotemp Flammable Material Undercounter Storage Refrigerator



5.5 cu. ft. (155.7 L) capacity; spark-free interior to prevent ignition risks

Rating: Not Rated Yet

Price

Sales price \$1,950.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

The Fisherbrand Isotemp Flammable Material Undercounter Storage Refrigerator (Model GTFBH05RFSA) is a compact, safety-focused solution for storing volatile materials in laboratory settings. With a 5.5 cu. ft. (155.7 L) capacity, it features a spark-free interior to prevent ignition risks, a temperature range of 1°C to 12°C, and a mechanical thermostat located externally for safe operation. Its enameled steel exterior, ABS plastic interior, and magnetic vinyl gaskets ensure durability and a secure seal, making it ideal for undercounter use in labs handling flammable substances.

Specifications:

Model: GTFBH05RFSA

Type: Flammable Material Storage (FMS) Refrigerator

Capacity: 5.5 cu. ft. (155.7 L)

Temperature Range: 1°C to 12°C

Refrigerant: R600a (Green, natural HC)

Defrost: Manual

Control Type: Mechanical thermostat (exterior-mounted)

Door Style: Hinged, 1 door (1 compartment)

Shelves: 4 adjustable

Exterior Finish: Enameled steel

Interior Material: ABS plastic

Gaskets: Magnetic vinyl for positive door seal

Evaporator: Aluminum coil

Compressor: Bottom-mounted

Dimensions (L x W x H) Exterior: 29 x 23.63 x 33.5 in. (73.6 x 60 x 85 cm)

Dimensions (D x W x H) Interior: 17.5 x 20 x 27.75 in. (44.5 x 50.8 x 70.5 cm)

Weight: 80 lb (36.2 kg)

Voltage: 115 V

Frequency: 60 Hz

Plug Type: NEMA 5-15 P

Additional Features: No internal electrical components, 4 in. clearance on sides/top, 3 in. at back

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com