

Fisherbrand Microplate Centrifuge Model # 14955300



850 RCF; 2500 RPM

Rating: Not Rated Yet

Price

Sales price \$750.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

For fast spindowns of PCR Microplates; small footprint with easy access removable rotor for easier cleaning. Easily spins down droplets and condensation in PCR Plates; easy insertion and removal of plates. Swing-out rotor design eliminates need for plate seals.

Microplate Centrifuge specifications

Capacity: (2) Standard PCR Microplates

Display: Digital

Max. Speed: 2,500 rpm

Max. RCF: 850 x g

Run Time: 1 sec. to 10 min.

Dimensions (W x D x H): 11.4" x 14" x 5.5"

Net Weight: 8 lbs

Electronic Brake

Electrical: 115 V, 50/60 Hz

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com