

## Fisherbrand Nutating Mixer - Variable Speed



## Shakers/Stirrers/Hotplates: Fisherbrand Nutating Mixer - Variable Speed

---

12 x 10.5 in. platform accommodates a variety of tubes and vessels

Rating: Not Rated Yet

### Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

### Description

Model: 88861043

#### Description

Variable speed from 2 to 60rpm offers gentle continuous mixing of samples

Low-profile, compact design enables placement at every workstation

12 x 10.5 in. platform accommodates a variety of tubes and vessels

Display will show last used settings, even after power has been turned off.

Operating temperature range - 5° to 40°C, usable in non-CO2 incubators, refrigerators and chambers

Continue to run by last set speed and the remaining time when the power recovered.

Continuous or timed operation with automatic switch-off

Audible and visual alarm - In timed mode, when the time reaches zero, alarm will sound five times, then the time window shows "End" and flickers all the time.

Alarm can be stopped by touch any key

Tilt angle: 20°

100-240V; 50/60 Hz; 1A

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)