

## Linea Schuster Foster Plus 22 L-N Autoclave





22L; 2.8" LCD display – 58x43mm.

Rating: Not Rated Yet

**Price**

Sales price \$2,500.00

[Ask a question about this product](#)

Description

6 sterilization programs

- Packaged Instrumental
- Unpacked Instrumental
- Plastic and Cotton
- Surgical Kit
- Fabrics
- Liquids

Total Cycle: \*60 minute;

\*Average time, subject to variation of  $\pm$  10%

Sterilization at 134° C and drying with the door closed. Packaged material.

Specifications

Autoclave Model: Foster Plus 22L-N

Autoclave dimensions: Depth: 62.0 cm x Width: 43.9 cm x Height: 47.5 cm

Chamber dimensions: Depth: 45.0 cm x Diameter: 24.9 cm

Chamber material: Stainless steel AISI 304

Capacity: 22 Liters

Sterilization: 121°C and 134°C

Rated voltage: AC110V-240V

Power: 1800 VA

Rated frequency: 50/60 Hz

Max pressure: 2.4 bar

Mfg date: 01.2020

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)