

Gast High Capacity Vacuum Pump Model DOA-P704-AA



1/8 hp, 115V AC, 25.5 in Hg Max Vacuum, 60 psi Max Continuous Pressure

Rating: Not Rated Yet

Price

Sales price \$250.00

[Ask a question about this product](#)

Manufacturer [Gast](#)

Description

This diaphragm combination compressor and vacuum unit runs continuously, without oil or excess noise, and provides light-duty vacuum or air pressure.

Specifications:

Manufacturer Part Number: DOA-P704-AA

Compressor Type: Diaphragm

Current Rating: 4.2 A
Frequency: 60 Hz
Inlet Connection Type: NPT
Inlet Size: 3/8 in
Maximum Continuous Pressure: 60 psi
Maximum Flow Rate: 1.9 cfm @ 0 psi
Maximum Intermittent Pressure: 60 psi
Maximum Operating Temperature: 104 °F
Maximum Speed: 1,575 RPM
Maximum Vacuum: 25.5 in Hg
Minimum Operating Temperature: 41 °F
Outlet Connection Type: NPT
Outlet Size: 3/8 in
Output Power: 1/8 hp
Overall Height: 8 in
Overall Length: 7.62 in
Overall Width: 5.12 in
Sound Level: 70 dBA
Voltage: 115V AC
Weight: 14.5 lb

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com