

Harvard Apparatus Model 11 Plus



Harvard Apparatus Model 11 Plus dual syringe pump, demo model

Rating: Not Rated Yet

Price

Sales price \$975.00

[Ask a question about this product](#)

Manufacturer [Harvard Apparatus](#)

Description

This Harvard Apparatus Model 11 Plus dual syringe infusion pump combines smoother flow and updated features to create a high performance syringe pump with a basic syringe pump price. New micro-stepping pump profiles deliver very smooth and consistent flow over the entire flow rate range. These pumps now feature non-volatile memory. The pump remembers its last syringe size, flow rate used and configuration settings. An advanced universal power supply means there is no need to change AC line switches, fuses or wires. Demo Model.