

Harvey MC10



Harvey MC10 benchtop microprocessor controlled autoclave, 10" x 15.5" chamber

Rating: Not Rated Yet

Price

Sales price \$2,850.00

[Ask a question about this product](#)

Manufacturer [Harvey](#)

Description

!!!SOLD!!!

This Harvey MC10 benchtop microprocessor controlled autoclave has a 10" x 15.5" chamber and permits one-touch operation from start to finish, and a time-temperature interlock that assures correct chamber conditions prior to cycle start-up. Drying begins automatically at the end of the sterilization cycle.

Please call us at (650) 871-5707 to check availability of this product!