

## Heidolph Basis Hei-VAP HL with Diagonal Coil Condenser



Temperature range (°C): 20 - 210; 4.5 L Capacity.

Rating: Not Rated Yet

**Price**

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Heidolph](#)

### Description

Rotary evaporator with Diagonal Coil Condenser.

Model: Heidolph Hei-VAP Value Digital

### Specifications:

Temperature range (°C): 20 - 210

Rotation speed: 10 - 280 1/min<sup>1</sup>

Heating capacity: 1400 W

Cooling surface: 1400 cm<sup>2</sup>

Bath capacity: 4.5 L

Lift distance: 155 mm

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)