

Heidolph PD 5201 Pump Drive



0.3 to 790mL/min flow rate with single-channel pump head

Rating: Not Rated Yet

Price

Sales price \$1,150.00

[Ask a question about this product](#)

Manufacturer [Heidolph](#)

Description

Easy to use and efficient for reproducible results and liquid transfer.

- Flow rates of 0.3 to 790mL/min with single-channel pump head
- Individual calibration of flow volume and flow rate
- Button for clockwise or counterclockwise pump operation
- Digital read-out of speed (rpm), tubing diameter, flow rate in mL/min., dosing volume, interval dosing, pause function
- Features electronic speed control and accuracy of $\pm 0.5\%$
- Low speed range from 5 to 120rpm
- Also suitable for multi-channel operation achieving flow rates of 0.005 to 320mL/min

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com