


Hirayama HV-110



Hirayama HV-110 toploading autoclave

Rating: Not Rated Yet

Price

Sales price \$7,650.00

[Ask a question about this product](#)

Manufacturer [Hirayama](#)

Description

!!!SOLD!!!

This Hirayama HV-110 toploading autoclave is a fully automatic, stand-alone, microprocessor-controlled, programmable system. It incorporates a built-in steam generator and plug into a standard 220V power outlet, thereby eliminating costly set-up charges. Water for steam generation is poured directly into the chamber. After sterilization, steam is vented automatically into a plastic bottle in front of the autoclave. This model has a 110L capacity and the cylindrical chamber is 16.5" diameter by 31.3" deep. Comes w/2 stackable baskets. Current model. Lists for over \$20K new.

Please call us at (650) 871-5707 to check availability of this product!