

## IEC Centra CL3 Centrifuge



IEC Centra CL3 Centrifuge (Rotor not included)

Rating: Not Rated Yet

**Price**

Sales price \$1,900.00

[Ask a question about this product](#)

Manufacturer [IEC](#)

### Description

The IEC Centra CL3 Centrifuge:

- Store and lock up to 99 protocols
- Speed and accuracies to the nearest 50 rpm as well as automatic speed calculation of the g-force
- Temperature with 1 °C increments
- Capable of reaching 4 °C from 20 °C in ~7 minutes with rapid condition means you always spin at the required temperature.

The CL3 series accommodates up to 1 liter of sample (4 x 250mL) accepting a wide variety of popular tubes including 15 mL, 50 mL, 225 mL conical, Vacutainer® tubes to 15 mL (16 x 125mm), microtubes, microplates (standard and deepwell), microscope slides, and an assortment of filtration devices.

### New Quick Release Rotor

The Thermo IEC 243 Swinging Bucket Rotor, utilizes convenient sealed Aerocarrier Adapters which accommodate sample volumes ranging from 0.25 mL to 250 mL.

**\*\*Centrifuge does not include rotor\*\***