

IKA Matrix Thermomixer Orbital Delta F1.5



Speeds of up to 1,500 rpm and shaker diameter of 3 mm.

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

Description

The IKA MATRIX Orbital ? F1.5 mixes and heats with a fixed insert for 24 x 1,5 ml tubes. Blood samples, pharmaceutical agents, DNA/RNA samples or ELISA assays – samples with even the smallest volumes -- mix reliably and completely.

Large and clear display with a convenient menu structure allows for intuitive operation; speeds of up to 1,500 rpm and shaker diameter of 3 mm ensure the liquid mixes quickly and effectively.

Programmable ramps allow for automatic execution of mixing and heating tasks; smart speed adjustment according to insert, attachment and filling volume. Aluminum die cast housing ensures secure positioning, excellent stability and everyday robustness.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax

sales@biotechequipmentsales.com