

## IKA MTS 2/4 Digital Four-Place Titerplate Shaker



4-place

Rating: Not Rated Yet

**Price**

Sales price \$400.00

[Ask a question about this product](#)

Manufacturer [IKA](#)

### Description

Features include: electronic speed control; digital timer; and alarm to indicate set time has expired.

### Specifications:

Type of Movement: Orbital  
Shaker Stroke: 3mm

Permissible shaking weight (incl. attachment): 1 kg  
Motor Rating Input: 35W  
Motor Rating Output: 13.5W  
Permissible ON Time: 100%  
Speed Min.: 50 rpm  
Speed Range: 0 – 1100 rpm  
Time Setting Range: 1 – 99 minutes  
Operating Mode: Timer and Continuous Operation  
Shaking Table Lock: No  
Working with Microtiter Plates: Yes

## Shakers/Stirrers/Hotplates: IKA MTS 2/4 Digital Four-Place Titerplate Shaker

---

Number of Microtiter Plates: 4  
Dimensions (W x H x D): 7.23? W x 12.59? D x 4.13? H (18.5 x 32 x 10.5 cm)  
Weight: 5.95 lbs. (2.7 kg)

Permissible ambient temperature: 5 - 40 °C  
Permissible relative humidity: 80%  
Voltage: 100 - 120 V  
Frequency: 50/60 Hz  
Power input: 45 W

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)